


DÉPARTEMENT DE PHYSIQUE  
DE L'ÉCOLE NORMALE SUPÉRIEURE  
LABORATOIRE KASTLER BROSSEL



# Manual for the computer control of LITHIUM

by Florian Schreck  
7.8.2001

# Contents

<b>1</b>	<b>Basic hardware knowledge</b>	<b>1</b>
1.1	Cards in the computer . . . . .	1
1.2	The principle of the bussystem Multiout . . . . .	2
<b>2</b>	<b>Lowlevel structures</b>	<b>3</b>
2.1	IO units . . . . .	3
2.2	Capsulated IO . . . . .	4
2.3	Renaming of ports . . . . .	5
2.4	IEEE communication . . . . .	6
2.5	Serial communication . . . . .	6
2.6	Timing . . . . .	8
2.7	The idea behind Portlist . . . . .	9
2.8	The idea behind Sequence . . . . .	9
2.9	Principle of an experimental run . . . . .	11
2.10	The RF sweeps . . . . .	12
2.11	Debugging . . . . .	13
<b>3</b>	<b>Mediumlevel structures</b>	<b>14</b>
3.1	The use of the Menu System . . . . .	14
3.2	The parameter administration . . . . .	16
3.2.1	Utilities for parameter administration . . . . .	18
<b>4</b>	<b>Highlevel structures</b>	<b>20</b>
4.1	The principle of an experimental serie . . . . .	20

4.2	The principle of a measurement queue . . . . .	21
4.3	Some special functions of the program . . . . .	21
4.4	Last words . . . . .	21
5	Basic modifications of the program	<b>23</b>
5.1	Adding a shutter . . . . .	23
5.2	Adding a box to multiout . . . . .	23
5.3	Adding a new Block in the experimental sequence . . . . .	23
A	The units of LITHIUM	<b>24</b>

# Chapter 1

## Basic hardware knowledge

### 1.1 Cards in the computer

Four cards take over special functions in the computer

- DAS1601

This card contains the  $100kHz$  timer chip used for the timing of the experimental sequence (a standard Intel clock chip, as also used for the internal clock of the computer). It contains a digital I/O chip with 24bit (also standard intel, as the one for the parallel port). Four digital inputs, two analog outputs and six analog inputs. The 24bit digital output is used for the bus system, one analog input for the Fluorescence measurement, one digital input for the Water and Temperature Error signal. The card is connected to a output box, made by our electronics. If another analog input has to be used, a OP driver has to be put into this box, the driver that are in place are too weak.

- DDA08

This card contains another 24bit of digital outputs and 6 analog outputs. If the DAS1601 breaks, the bus can be hooked up here, only the adresses in Multiout.pas have to be changed. These outputs are thought as spare outputs in case no output is left on the bus and there is no time to installe a new output box. But you should better use the bus system.

- IEEE card

This is for communication with IEEE gear, like synthesizers, . Useful tools came with this card and are in `c:\debug488\` .

- serial card

Two serial ports for the communication with the Nd:YAG and the data acquisition computer. Another two ports can easily be added by adding another card (already tried), but pay attention to the interrupts and IO addresses.

## 1.2 The principle of the bussystem Multiout

The bus system (Multiout) has a 16 bit data bus, 6bit adress bus and a 1bit subbusselection and a 1bit strobe bit. The bus is hooked on a driver box, splitting it up on two cables (the subbusses). At the moment always both subbusses are activated and it is rather stupid to sparate them, I would change the electronics of the driver box to use the subbusadressbit as a normal adress bit. Boxes with analog or digital outputs are hooked on the bus. Each box is optoisolated. Each 16bit latch or 16bit DA converter has its own adress, decoded by a 8bit comparator. Most of the boxes decode 8adress bits (one of them is strobe) the oldest boxes only five, they have to be modified if new boxes are hooked up since the 32 5bit addresses are used up. The modification is very simple and consist of wirering the 3 highest bits from the output of the optocouplers to the input of the adress decoder. The adress is coded by hooking the 8bits of the second input of the comparator to GND or 5Volts. Pay attention, the bits are rather mixed up (21436587 in stead of 12345678). Some of the boxes have an activated error input. The outputs can only be programmed if a 5V signal enters there. If it goes to GND the outputs go to GND too. The digital output box for the IGBTs and the Cameras is split in two, 8bits are secured, 8bits not. This is to prevent that all IGBTs go off at the same time which can provoque overtension.

# Chapter 2

## Lowlevel structures

Now I will describe the ideas behind LITHIUM and the relationships between the different units.

### 2.1 IO units

Four units, for DAS1601, DDA06, Multiout and the printer port LPT speak with the in/out ports, each contains an initialisation routine and input output routines, see exapmle.

```
unit DAS1601;
unit DDA06;
unit Multiout;
unit LPT;

PROCEDURE DAS1601Init;
procedure DAS1601analogOut(Channel: byte; Value: word );
{channel:0..1; value:0..4095}
Function DAS1601analogIn(channel:byte):word;
{channel:0..7;value:0..4095}
Procedure DAS1601digitalOut(channel:byte;Value:byte);
{channel:0..3} Function DAS1601digitalIn(bitNR:byte):boolean;
{bitNr: 0..3} Function DAS1601FastanalogIn:word;
```

In addition these units contain code necessary for the DEBUG option, writing all the IO traffic to a file on the harddisk if necessary. A speciality

of MULTIOUT.PAS is that here the channel numbers are translated to the physical addresses, using these tables (look them up in MULTIOUT.PAS for completeness :

```

                                {0, 1, 2}
Const DigAdr:array[0..NumbDig-1] of byte=(29,30,31);
    AnaAdr:array[0..NumbAna-1] of byte=(20,21,22,23,24,25,...
                                {0, 1, 2, 3, 4, 5,...

```

The numbers in { } give the addresses from the program calling MULTIOUT.PAS point of view, the numbers in ( ) give the physical address. If a new box is added these tables have to be extended. (First increase NumbDig or NumbAna, then add unique addresses. Each analog box has four DACs and so four addresses.)

## 2.2 Capsulated IO

From the programmers point of view, it would be nice to have one addressspace for all cards, so he does not always have to think to which card exactly he has to write when a bit is set. That's one of the functions of IO.PAS. The other is that several write commands going to the same digital 16bit port, but to different bits can be accumulated. The commands are only stored in a buffer containing the value of the port it should have. Only when a IOUpdate command is found, IO searches through its list of ports and writes the ports which were changed. This increases speed and makes it possible to change several bits on one port exactly at the same moment. That's used for the switchoff of the Magnetic trap and the triggering of the camera system for example. A WAIT command always calls IOUpdate before waiting. The table linking the virtual address space to the boxes address space (for Multiout that's NOT the physical address space, see before) has the following structure (given is only the digital ports table, the analog ports one is more or less the same, see IO.PAS):

```

const DigBuf:array[0..MaxDigOut-1] of TDigBuf=(
    (card:CMultiOut;channel:0;NrBits:16;value:0;update:true), {channel 0}
    (card:CMultiOut;channel:1;NrBits:16;value:0;update:true), {channel 1}

```

```

(card:CMultiOut;channel:2;NrBits:16;value:0;update:true), {channel 2}
(card:CDDA06 ;channel:0;NrBits:8;value:0;update:true), {channel 3}
(card:CDDA06 ;channel:1;NrBits:8;value:0;update:true), {channel 4}
(card:CDDA06 ;channel:2;NrBits:8;value:0;update:true), {channel 5}
(card:CDAS1601 ;channel:3;NrBits:4;value:0;update:true), {channel 6}
(card:CLPT ;channel:1;NrBits:8;value:0;update:true) {channel 7}
);

```

On the righthandside the addresses from the rest of the programs point of view are given. Left the outputsystem and the channel of this system is given. Space is reserved for the buffer and a statusflag is set. NrBits contains how many bits a channel has. This table has also to be extended if a new box is added. The rest of the program (which means EXPIO.PAS) speaks with the ports only through IO.PAS. In addition IO.PAS contains some code necessary for a correct handling of Portlists (see there), storing and recalling the buffers.

## 2.3 Renaming of ports

It is incomprehensible if in the final experimental sequences code the addresses of special bits appear. If for convenience an other outputbit should be used for a function, the whole code would have to be searched for the usage of that bit. That's why every port gets its own name describing as clearly as possible its function, also if the name gets very long (copy paste helps). Here is an example:

```

procedure ShutterProbe(onoff:boolean);
begin
  if SequenceOut then IOdigitalOut(1,9,not onoff)
  else SequenceBooleanArg(ShutterProbe,onoff,0);
end;

```

A call of ShutterProbe calls bit 9 of address 1 (pay attention, bits here are numbered from 0..15, on the boxes they are numbered from 1..16.) The experimental sequence is executed twice for each experimental run. The first time to gather general information (like total time consumption), then



for real. That's why the outputs are blocked if SequenceOut is false. In fact then the command itself is stored in a list in memory and is called when this list is worked through. The list itself serves to displace the actual output routine from their place when it is called during the sequence (eg. give the command to a shutter earlier). See Sequence.pas for more info.

Other things done by ExpIO.pas are implementing simple IEEE commands, organising the execution of WAIT commands for Portlist and Sequence (see there).

## 2.4 IEEE communication

The unit IEEE contains routines for this. Each IEEE machine has its own supplementary unit, like NDYAG.PAS, SIGGEN.PAS, PULSGEN.PAS hiding the IEEE commands behind nice names. SIGGEN.PAS is a bit special because it implements an abstract synthesizer which is mapped to the actually used one, so that one synthesizer can be exchanged with another just by changing the type in the list of synthesizers in SIGGEN.PAS. The supported types are: SRS, Marconi, Rhode&Schwartz and the old Synthesizer from Grynbergs lab. The GPIB addresses of Synthesizers are given in SIGGEN.PAS, all the others in EXPIO.PAS.

## 2.5 Serial communication

Two units are at disposition for serial communications: serial1.pas and serial2.pas. One uses DOS interrupts and is slow, but works. The other speaks directly to the chip is fast, but when I tried it, there were many missing characters after transmission. So an error correction protocol would have to be implemented to use this unit. Now the simple DOS unit is used for NDYAG.PAS and VISION.PAS. NDYAG.PAS implements all possible commands for the NDYAG. The COM port used has to be set before calling this unit using NdYAGSetComPort(COMx) (x=1,2). This is done during the initialisation in EXP.PAS.

The communication with Vision is done through VISION.PAS. The com port is set there. If new parameters which Vision has to know before the image is taken, like camera parameters, have to be transmitted this unit

has to be modified. But pay attention, if the communication grows longer, the VisionWaitingTimes have to be increased in the configuration menu. The main commands in VISION.PAS are

```
procedure VisionSetCameraParameters(...
function VisionCameraReady(Camera:byte):boolean;
function VisionStartSerie(...
function VisionStopSerie:boolean;
function VisionTakeAbsorptionPicture(...
procedure VisionSendPictureData(
procedure VisionMessage(Message:string);
```

The communication goes as follows. First the Camera is initialised with SetCameraParameters. When CameraReady responds with true, eventually a Serie is started. The experimental sequence is executed a first time to determine the total time. This together with some other parameters is sent with VisionTakeAbsorptionPicture. The returned # READY\* puts both computers in synchronisation. Fine tuning is done with the AndorTriggerCalibrationAbsorption and AndorTriggerCalibrationAbsorption parameters accessible in the Andor menu. For the Matrox camera this is not necessary. 100ms before the image is taken, Vision puts the Andor system in the external triggerable mode. In this mode the CCD isn't cleaned well any more, that's why it's invoked as late as possible and that's why the total time has to be known before the experimental sequence starts. VISION and LITHIUM synchronise each other with some Cameratriggers and # READY\* signals sent back to take the three images (or one if the noise images comes from harddisk). After finishing of the image taking and digitising VisionSendPictureData sends supplementary information and Parameters modified since the last image. With VisionMessage messages can be sent to Vision which are shown in Visions statusline and in the statusline of all Lifeconnected computers. This serves to see at a glance where in a sequence the experiment is, when sitting in the office and to see if a Temperature or water error did appear. Some additional commands can guide vision to do any of the commands also available in Visions Menus.

## 2.6 Timing

TIMING.PAS implements a good delay routine using the DAS1601 Counter 0. If the DAS1601 card isn't detected the PC internal timer is used. The timer is a 16bit timer counting at a  $100kHz$  rate. Initially upon call of StartCounter, its value is set to \$FFFF and also stored as last value in the OldCounter variable. Then it starts to count down. When a WAIT command is met, first IOUpdate is invoked to write all outputs. This and all the operations since the last WAIT command take an unknown amount of time. Now the new value expected for the counter at the end of the wait is determined by subtracting the number of counter cycles corresponding to the waitingtime from OldCounter. If the counter is already below this value, the program took too long. Synchronisation is no longer guaranteed and an error printed on the screen, indicating the identity of the too short wait. This musn't be catastrophic, in most cases the sequence is still fine, but a bit longer. If the counter is higher, TIMING waits till it reaches the value calculated, stores this value in OldCounter and gives control back to the calling routine. Like this some IO operations can be added without influencing the timing at all, the sequence becomes completely repeatable (at a timescale of  $10\mu s$ ). The most difficult thing to take into account has been the timeroverflow appearing each  $500ms$ . This means that a Wait command has to appear at least every  $500ms$ . (This can be circumvented by calling CheckCounter from time to time, which calculated the time elapsed since the last wait and stores the actual countervalue. This is taken into account when the next wait appears.)

You have to remember to call startCounter shortly (less than the first wait-time and less than half a second) before the first Wait command. This is necessary if Waits are used out the experimental sequence, like in the menu. To identify Wait commands in the sequence to each delay time a femtosecond delay is added which of course doesn't play a role. But the amount of femtoseconds identifies the wait. Example:

```
WAIT(100+ID*42);
```

The function GetTimeSinceLastWait:double; is used when fine synchronising the two computers and gives in that case the time needed by the Vision system to get into the gears.

## 2.7 The idea behind Portlist

Sometimes commands to the ports have to be emitted in a very fast sequence. When using the normal way of writing using IO.PAS and afterwards IOUpdate, about  $70\mu s$  of administration overhead are spent each IOUpdate. When using a Portlist, the commands are not executed at once. Rather the values really written to the IO ports of the computer and the Waits between them are stored in a list. During the creation of this list the administration time is spent. Next the list is put to the ports without "thinking". In the following an example taken from MAGTRAP.PAS, the capture in the Ioffe trap:

```
SequenceStartistDuringWait(159*ID+2);
SwitchMOTDecharging(off);
IOUpdate;
SetMOTVoltage(PinchCapacitorVoltageStart);
SetMOTCurrent(MaxMOTCurrent);
Wait(160*ID+PinchCapacitorChargingTime);
SwitchIGBT(1,off);
SwitchIGBT(2,off);
SwitchIGBT(3,off);
SwitchIGBT(4,on);
.
.
.
SequenceStopPortList;
```

Its the  $159ms$  wait, that is used to do all the administrative work. After its done the rest of the  $159ms$  is waited and then the Portlist is executed. Why the commands are preceded by "Sequence" is explained in the next section. This is not necessary under all conditions.

## 2.8 The idea behind Sequence

The initial idea behind sequence was the following. A shutter has to be opened. But when at the point in the sequence when you need it to open,

it's too late to give the command since the shutter has a delay. You have to write to the output port BEFORE the command appears in the sequence. This is possible when first storing the whole experimental sequence in memory. If a command appears that has to be executed earlier, its position in the list is just changed. This works, but there are many problems in the details. The memory of our computer is not sufficient to store every single command. So you can regroup commands in subroutines and only store the call to this subroutine. This is essential when doing ramps of currents or radiofrequencies, when thousands of commands are executed. Inside these blocks regrouping is of course impossible. If a command appears just after a block which should be preexecuted just a few milliseconds before, that's not possible. Thus a sufficiently long security wait should be put in place before any command using pretiming is used. Next, to make work portlist and pretiming work together was a big effort. In fact in daily life it's so inconvenient to use pretiming that I've abandoned it. The problem is before all that it's extremely difficult to debug an experimental sequence when the commands can hop around.

But sequencing has another advantage. That's that the experimental sequence is executed twice. Like that one can collect global parameters of the sequence like the "life of the atoms" (a descriptive string marking with two characters each step taken by the atoms) or the total time. These parameters can be transferred to VISION before experiment is started. For the total time that's necessary. Of course these informations can also be obtained without storage of the commands in memory. So in principle one could modify the program in a way that the whole experimental sequence is stored as one block. This will help if one day the memory becomes short. An other point where the preexecution is used is the programming of the puls generator. Its programmed during the preexecution and its different setups (up to 12) are stored in its memory. During the experimental sequence these settings are just recalled, which is much faster.

Now some examples of the use of the sequencing commands:

```
procedure ShutterProbe(onoff:boolean); begin
  if SequenceOut then IOdigitalOut(1,9,not onoff)
  else SequenceBooleanArg(ShutterProbe,onoff,0);
end;
```

The command `SequenceBooleanArg` stores `ShutterProbe` on the heap, together with its argument `onoff`. The delay is 0 and would be negative for a preexecution and positive for a later execution than the place at which `ShutterProbe` is encountered in the Sequence. The procedure `ShutterProbe` has to be called "far", which is automatically the case when it is declared in the interface of a unit like `EXPIO.PAS`. In other cases the `far` modifier has to be put after the declaration of the procedure. Since PASCAL can't analyse at runtime the arguments of a procedure to which it gets a pointer reference, different storage commands had to be implemented for each possible argument list (see `SEQUENCE.PAS`).

## 2.9 Principle of an experimental run

Now let's take a look at an experimental run, paying attention how Sequence is used. I think the code should by now be self explaining. You can find it in `EXP.PAS`, in the routine `TakeAbsPicture`. Probably its a good exercise to go through this routine and the whole experimental sequence to understand how it is programmed by using all the lower level structures described before. That's how the experiment is programmed in the daily life. By looking how things are made, copying and modifying them you should be able to do most of the things necessary for the future of the experiment.

```
StartCounter;
ReInitialiseSystem;
StartSequenceDuringWait(62*ID+334);
  SynchronizeToLine; syncs whole sequence to line trigger
  do stuff which will be stored on heap
  for example the whole experimental sequence
  if DoTransferQuadrupoleExpansion
  then TransferQuadrupoleExpansion; just a example
  cut trap for TOF
  PictureWait:=GetTotalWait; get time of the sequence till picture
  trigger camera
PauseSequence;
  things are executed at once now.
  Check Fluorescence
```

```

UnPauseSequence;
StartCounter;
  VisionTakeAbsorptionPicture(PictureWait-AndorTriggerDelay-
    GetTimeSinceLastWait-AndorTriggerCalibrationAbsorption,
    Atompath,...
StopSequence;  and off it goes! the stored sequence is executed
Cleanup: ramp down power supplies, sent additional parameters to
Vision VisionSendPictureData(...)

```

And this is how a sequence block is implemented. There are lots of them in MAGTRAP.PAS. Ramp(command, StartValue, EndValue, Time, CommandTime) is a very useful command to do ramps. The command is executed repeatedly with arguments going linearly from StartValue to EndValue. The whole ramp takes Time. Every CommandTime, command is called. CommandTime is typically 0.2, but can be larger if command regroups several commands writing on the ports. AtomPath is the string describing the sequence.

```

procedure TransferQuadrupolExpansion;far;
begin
  FinalTrap:='transferred MagTrap with Quadrupol Expansion';
  AtomPath:=AtomPath+'Qe';
  SequenceStartBlockNoArg(TransferQuadrupolExpansion);
  Ramp(SetTransferCurrentControlVoltage,TransferCurrentEnd,
    TransferQuadrupolEndCurrent,
    TransferQuadrupolExpansionTime,0.2);
  SequenceEndBlock;
end;

```

## 2.10 The RF sweeps

The unit RFSWEEP.PAS contains a class describing RF sweeps. It contains methods to save and load a sweep, to display it graphically, to create a Menu to modify it by hand and to send it to the Synthesizers.

## 2.11 Debugging

To debug the experimental sequence, all operations occurring during the sequence can be written to a protocol file. There are different modi for this. One creates a table that can be displayed with Origin, giving a multichannel oscilloscopes view of the experiment. Another gives a ASCII list of all output operations. DEBUG.PAS contains some Basic FileI/O functions for this. The routines really writing in the file are in the DDA08.PAS, Multiout.pas, the SERIAL.PAS and the IEEE.PAS units.



# Chapter 3

## Mediumlevel structures

Now we come to the structures concerning not the experimental sequence but the user interface and the parameter administration.

### 3.1 The use of the Menu System

MENU.PAS contains procedures to declare either a Menu, giving the option to start subroutines, which can also describe submenus, giving so rise to a menu structure, or Dialogs, which are used to modify the experimental parameters. When the user is doing nothing IdleProc is called which is mapped to an Idle routine in EXP.PAS showing the clock, beeping at 7:50 pm for the pot des elevés, closing the atomic beam shutter if the experiment is not used for quite a while, showing the error status and starting the screensaver. I won't explain MENU.PAS. Its application becomes clear from looking at our menustructure implemented in MENUSYS.PAS. Here an example for a Menu, calling utilities from EXP.PAS. If this menu calls an equal procedure describing another menu the user get the impression of a menu structure. The Menus which have to be adapted when the experiment grows are the manual operation menus. Here actions are executed when a button is pressed. Don't forget StartCounter at the beginning of such a small experimental sequence if WAITs are used and an IOUpdate at the end.

```
procedure MenuUtils1;  
var MenuInfo    : pMenuInfo;
```

```

    iMain,iTest7LiExpInjection, and so on
    iSweepDigitalOut,iInit,iCalFluo,Choice:integer;
begin
  NewMenu(MenuInfo);
  AddChoice(yellow,' ',yellow,MenuInfo,'Initialise System [I]',
    'Init System','I',iInit);
  AddChoice(yellow,' ',yellow,MenuInfo,'Test 7Li Injection [J]',
    'Test 7Li Injection','J',iTest7LiExpInjection);
and so on
  AddChoice(yellow,' ',yellow,MenuInfo,'Main Menu [M]',
    'Main Menu','M',iMain);
  choice:=-1;
  repeat
    DoMenu(MenuInfo,' Experiment Utility Menu ',Choice);
    if Choice=iProbe then TakeProbePicture
    else if Choice=iInit then InitExp(true)
    end else if Choice=iprobeoff then begin
      SwitchProbeBeamGPIB(off);
      ExposureTrigger(off);
      IOUpdate;
    end else ...
and so on
  until (Choice=0) or (Choice=iMain);
  DeleteMenu(MenuInfo);
end;

```

And here the example of a Dialog modifying the double parameter Frequency defined elsewhere. Normally its unnecessary to define dialogs by hand since the PARAMS.PAS unit in collaboration with the REGISTER.PAS unit does this automatically.

```

procedure BlinkParams;
var InfoDialog : PDialogInfo;
    UpdateDialog : boolean;
begin
  NewDialog(InfoDialog);

```

```

AddT(yellow,InfoDialog,'Blink Parameters');
AddR(yellow,InfoDialog,Frequency,'Frequency [Hz]',
      'Frequency',0,10000,1);
BackupVars(InfoDialog);
DoDialog(InfoDialog,'Blink Parameters',UpdateDialog);
end;

```

## 3.2 The parameter administration

An experimental sequence is controlled by a lot of parameters (at the moment about 700). All these parameters are loaded from disk when LITHIUM is started, stored to disk when it is quit or a serie of experimental runs is started. They are also sent to Vision and stored with each image. (In fact only the first time when Vision is started ALL the parameters are sent. Afterwards only the modified ones are sent, to gain time.) Clearly a system is necessary to handle all these parameters. The necessary routines for this system are in REGISTER.PAS (helped by STRING.PAS). Besides doing the things described before, the necessary Menuentries and Dialogwindows are automatically created. The actual declaration of the parameters takes place in the unit PARAMS.PAS. Each parameter has to be declared once in the interface part of the unit. Next it has to be registered in the implementation part. Here is an example:

```

unit params.pas
interface

const MaxCompressIoffeTrap=17;

var DoCompressIoffeTrap:array[1..MaxCompressIoffeTrap] of boolean;
    CompressIoffeTrapTime:array [1..MaxCompressIoffeTrap] of double;
    BarCurrentIoffeTrap:array [1..MaxCompressIoffeTrap] of double;

implementation
uses Register;

procedure RegisterParameters;

```

```

begin
  some code here
  DummyReg:=new(PRegistry,init(@Registry,Registry,true,
    'Doppler cooling Parameters',Dummy));

  RegisterBoolean(MarkColor,@DoCompressIoffeTrap[1], True ,
    'C1DoCompressIoffeTrap','Do Compress Ioffe Trap 1 ? [Yes/No]');
  RegisterDouble (StdColor,@CompressIoffeTrapTime[1],0,0,50000,
    'C1CompressIoffeTrapTime','Compress Ioffe Trap Time 1 [ms]');
  RegisterDouble (StdColor,@BarCurrentIoffeTrap[1],
    0 ,0,1000 , 'C1BarCurrentIoffeTrap',
    'Bar current in compressed Ioffe Trap 1 [0..1000 A]');
  some code here
end;

```

The first line in the RegisterParameters unit declares a new Parameter submenu. Registry is a pointer to the list of parameterpages of type PRegister. The boolean parameter tells REGISTER to show this parameters in the menu with the title 'Doppler cooling Parameters'. Dummy is a procedure without arguments, called whenever the corresponding parametermenu is quit. This is used for example on the MOT parameter page to update the AOM frequencies when the MOT detunings were changed.

Next the command RegisterBoolean registers a boolean variable deciding if the magnetic trap will be compressed. Since many deformations of the magnetic trap appear in the sequence, this is in fact the first boolean in an array of booleans declared in the implementation section. This is a very typical manner to simplify the experimental sequence. The procedure compressing the trap is called with a byte parameter, defining which of the 17 compressions should be performed. It is also typical to call a certain block of the experimental sequence in dependence of a DoItOrNot boolean variable. Lets look in detail to the RegisterBoolean routine. The first parameter is the color in which the dialogentry will appear. The color MarkColor is reserved for parameters deciding about the fate atoms will take. If a boolean parameter with MarkColor is switched on in a parametermenu, in the main menu a green spot will be displayed at the side of the call to the parameterssubmenu. This is very usefull to see at a glance

which options are switched on or off. The next parameter is a pointer to the boolean variable. Then its initial value follows, which is by the way rather unimportant, we could really get rid of that. For double parameters the parameter range would follow now. Next the unique name under which the parameter will be stored. Since only the first 20 or so characters are significant it is a good habit to number the names at the beginning if several of them exist in an array. The last parameter is the question asked the user in the dialog menu. That's it. In the next lines the other parameters for the compression of the trap are registered.

PARAMS.PAS contains also some code to store the RFSweeps.

A last small thing to know is that the SHUTDOWN.PAS unit saves a shutdown confirmation (SHUTDOWN.LOG) upon successful writing of the parameters when leaving the program. It is erased again when the program is loaded. If the program is quit anomalously it is not found when the program is started. This indicates that LITHIUM and VISION are not synchronised any more and all the parameters are transmitted again.

### 3.2.1 Utilities for parameter administration

As you know, its "chiant" to have to change a lot of parameters at a time. For this the SETUPS.PAS and CONFIG.PAS units were created. SETUPS.PAS reads the file SETUP.DAT in which subsets of all parameters are declared. A menu is shown with the names of these subsets. If this is selected the subset of parameters is loaded from the file. An example of the SETUP.DAT structure:

```
***CONFIGURATION***
```

```
Small Angle Optical System Calibration
```

```
AndorEichY=0.00306 ; Y Calibration [mm/Pixel]
```

```
AndorEichX=0.00306 ; X Calibration [mm/Pixel]
```

```
***CONFIGURATION***
```

```
Wide Angle Optical System Calibration
```

```
AndorEichY=0.0307 ; Y Calibration [mm/Pixel]
```

```
AndorEichX=0.0307 ; X Calibration [mm/Pixel]
```

**\*\*\*CONFIGURATION\*\*\*** is a keyword, declaring a new subset. The next line has to contain its name, then the parameterdeclarations follow. The example shows the use of this to change the calibration quickly between the two lens systems of the ANDOR camera.

The CONFIG.PAS utility stores and recalls all parameters. It resembles the save and loadgame option of a video game. About 20 slots are open. When selected in the savemenu, the user is asked to give the configuration a name and then all parameters are stored in this slot. Physically this corresponds to writing all setup files (RFSweep, Unimess and Parameters) into a subdirectory named SLOTA to SLOTZ. In addition a info.dat file is saved containing the name and creation date and time of the saved configuration. The contents of all these files is loaded when the Configmenu is build up, that's why it takes a small amount of time. Loading works analogous.

# Chapter 4

## Highlevel structures

The last units not mentioned until now will be described here. Lets take a top down look on the program. The main program is LITHIUM.PAS. It calls InitExp from EXP.PAS and MainMenu from MENUSYS.PAS. The experimental sequence is defined in EXP.PAS and MAGTRAP.PAS. In EXP.PAS are many more features, like the screensaver, experiment initialisation, all the test functions for analog and digital outputs and the serial ports test of laser injection, commands for the measurment queue, Fluorescence oscilloscope...

### 4.1 The principle of an experimental serie

The next feature to be described are experimental series. The procedures involved in this are declared in UNIMESS.PAS, EXP.PAS and MIX-DATA.PAS. UNIMESS.PAS administrates all the data necessary to describe a serie. It makes the menus necessary to create one, stores and loads their description in UNIMESS.DAT and finally calls DoUniMess in EXP.PAS. An universal measurement consists of executing the same experimental sequence over and over again by linearly varying up to four parameters. As you all use the program I won't describe more what it does. Since the measurements are executed with an arbitrary order, MIX-DATA.PAS provides a way of doing this in a quasi-random way. If the serie is stuck after having gone half through, it can be resumed, the starting point has to be changed changing MeasurementStartPoint in the Configuration parameters menu. Vision reorders the measurments again in the correct

order. A Vision command can be executed at the end of a Serie. By the way, the BigMeasurement option is in the meanwhile disfunctional.

## 4.2 The principle of a measurement queue

A serie only varies some parameter linearly. A measurement Queue is much more versatile. When creating a queue all parameters are stored on harddisk. An entry in the queue can consist of a simple experimental run or of a whole serie (not involving RF sweep parameters, which are treated a bit apart from the other parameters). When programing a queue it is useful to first test and then store the different parameter sets using the CONFIG utility. Next the old Queue is erased and then one by one the configurations recalled and stored with QueueAbsorptionImage. At the end the Queue can be executed once, or in circle. This makes it possible to test the stability of the system at two different places in the sequence. Or different Series can be started automatically one after the other. Physically Queues are stored in a special directory on the harddisk. Every element has its own subdirectory. All parameters and a file MESS.DAT describing the element are stored. It would be very easy to write an utility displaying the contents of the queue and making giving the possibility to selectively delete elements from the queue and not only to erase all.

## 4.3 Some special functions of the program

There are a lot of test functions for all kinds of in and outputs, the serial port, the IEEE port and lot of more stuff defined in EXP.PAS and accessible via the Utility menus. Just take a look and explore.

## 4.4 Last words

Some words to finish. Don't try to keep in mind the position of every routine in MAGTRAP.PAS or EXP.PAS. If you want to find something, first look it up in the parameter Menu of the program. You will see the name of the parameter used in Lithium in the status bar. Next search for



this parameter in EXP.PAS or MAGTRAP.PAS and you will find the point where you want to modify the program.

The Turbo Pascal editor has many very usefull, but somewhat hidden functions. You can find the in the Menu help-Content-editor. My favorite are Ctrl K 1 and Ctrl Q 1 to store the place of the cursor and hop back later (1 can be any number). Ctrl Y erases a whole line. Shift Cursor hops the cursor from word to word.

# Chapter 5

## Basic modifications of the program

### 5.1 Adding a shutter

In EXPIO.PAS a new routine has to be declared, copy a old shutter routine. Give it a new name. Don't forget to modify the output port and to adapt the name in the SequenceRegister command. Finally initialise the Output in Init in EXP.PAS or in the case of a shutter in PrepareShutters. A Shutter command has also to be copied to OpenShutters and CloseShutters. Finally the element has to be reinitialised in the procedure ReInitialise after a sequence.

### 5.2 Adding a box to multiout

Extend the tables in MUTIOUT.PAS and IO.PAS.

### 5.3 Adding a new Block in the experimental sequence

Look up the sections PARAMS.PAS and SEQUENCE.PAS in this manual.

# Appendix A

## The units of LITHIUM

program codetable;

Writes ASCII codes and characters on screen

unit config;

uses params,menu,dos,crt;

Stores and Loads configurations to slots. Configuration consists of Parameters declared in Params.pas Creates Menus necessary for these operations.

UNIT DAS1601;

Basic analog and digital In/Out routines for the DAS1601 card. Our bus system is connected to this card. The DAS1601 Timer is used in Timing.pas. The Fluorescence and Error signals are connected to this card too.

unit DDA06;

uses Portlist;

Basic analog and digital In/Out routines for the DAS1601 card. This is a spare card with 6 analog outputs and 24 digital outputs independent of the Bus system. In case the DAS1601 card breaks, the Bus can be hooked to these 24bits.

unit Debug;

Basic FileI/O functions for the debugging option. Debugging means writing all experiment outputs also in a file. This emulates a multi channel oscilloscope.

```
unit Exp;
```

```
uses unimess,IO,ExpIO,Params,Timing,Vision,MixData,DAS1601,  
      DDA06,MultiOut,Serial,IEEE,menu,menusys,MagTrap,Debug,  
      SigGen,Sequence,NdYAG;
```

This is the main program. It uses the different units to establish the connection with Vision, prepare an experimental serie or a simple experimental run. When the experimental sequence goes to the magnetic trap, MagTrap.pas is called, simply to keep the filesize of Exp.pas short. MagTrap.pas can be considered as part of Exp.pas. There are also other functions here: The screensaver, experiment initialisation, all the test functions for analog and digital outputs and the serial ports test of laser injection, commands for the measurment queue, Fluorescence oscilloscope, measurement queue

```
unit ExpIO;
```

```
uses Params,crt,IO,IEEE,Timing,Sequence,PortList,debug,Pulsgen;
```

This unit is used to give names which are easy to memorize to the different output and input ports. It also blocks the output if the experimental sequence is executet in the preparation mode and stores the command in the sequence list. It contains also simple IEEE Commands and all but the synthesizers GPIB adresses. The latter are in siggen.pas. It also reimplements some comands of Sequence.pas and contains a general Wait command that is reinpertreted in dependence of the actual execution state of the experimental sequence (preparation, portlist, SequenceOut). This is used also to calculate the total execution time.

```
unit HP8563A;
```

```
uses IEEE;
```

This unit isn't used in Lithium but can be used to read traces from the HP8563A READHP8563A is the mainprogram using this unit.

UNIT ieee;

uses debug;

This unit contains everything necessary to communicate using the IEEE interface. The communication can be written in a protocol in debug mode. Don't forget about the very useful

c:\DEBUG488\DEBUG488 (and c:\DEBUG488\IEEE488) utility, which monitors IEEE commands in real time, if you have trouble with IEEE.

unit IO;

uses crt,DAS1601,DDA06,MultiOut,LPT;

This unit unites all Input and output operations The digital ports are projected to virtual 16 Bit ports All channels are sorted and get a unique channel number The number of IO operations is minimized using buffers Data is written out when IOUpdate is invoked or if IOdigital/analogIn/OutDirect is used. If the calculation effort (30s per IOUpdate) is too long and slows down the experiment, a Portlist has to be used. A Portlist is constructed by executing the sequence without actually writing to the ports. All I/O operations and WAIT commands are stored in the list. Then the list is written out "without thinking", which is faster. The analog Values go from 0 to \ \$FFFF even for 12bit DACs. Look up the tables in this unit to get the link between channel numbers uses in expio.pas and channelnumbers used in the individual output system driver units DAS1601, DDA6 LPT and MultiOut.

program keypress;

helper program that translates keystrokes in keycodes.

UNIT LeCroy;

Utilities to read waveforms from a LeCroy 9314 or 9310 digital oscilloscope using its GPIB interface and a CEC GPIB card on the PC side. Not used in Lithium.pas.

```
program Lithium;
```

```
uses menusys,exp,expio,crt,io;
```

Mainprogram. Calls the experimental initialisation from exp.pas and then starts the Menusystem in Menusys.pas. Good place to put small test routines for output ports that one wants to use only once and deletes at once afterwards.

```
unit LPT;
```

Also the parallel port contains digital outputs. Once the DAS1601 which delivers the 24 bit for the bus system broke and we hooked the bus system to several parallel ports and the DDA06, which works perfectly well.

```
unit MagTrap;
```

```
uses IO,Params,Sequence,ExpIO,Exp,crt,timing,RFSweep,SigGen,Pulsgen;
```

Exp.pas calls the continuation of the experimental sequence here as soon as the atoms get optically pumped. Here are all the routines for optical pumping, capture, the transfer and capture in the Ioffe trap, Doppler cooling, calls to evaporation (RFSweep.pas), everything for internal state transfers, optical trap and Feshbach field.

```
UNIT Menu;
```

```
uses strings,CRT,dos;
```

The routines necessary to create Menus and Dialogs are implemented here. If the program is idle and waiting for user action it can call an idle routine if defined. In LITHIUM the Idle routine is in Exp.pas and tries to protect the experiment (closing atom beam shutter), shows clock, status, and starts the screensaver.

```
unit Menusys;
```

```
uses Register,crt,Menu,Params,Exp,IO,ExpIo,Vision,MagTrap,UniMess,Debug,  
MultiOut,Setups,strings,config,NdYAG;
```

This unit describes the menu system. An exception are automatically generated menus like for example the parameter list. Some special menus are the Manual operation menus. They contain

code to switch click-clacks IGBTs, program power supplies etc.

```
unit MixData;
```

a helper class to pseudorandomise the order of the parameters in a measurement series.

```
unit MultiOut;
```

```
uses PortList,Debug,crt;
```

Software interface to our homebrewn digital/analog output system. Look up the tables here to get the link between hardware addresses and the addresses uses in IO.PAS to access Multiout.pas. contains also code for debugging. Our bus contains 16 data bits, 7 adress bits and a strobe bit. To write to a box, data and adress are put to the bus. 10s afterwards when data and adress are stable the strobebit is flashed for 10s. This prevents accidental writing of data during the change of adress or data bits. Boxes with floating or nonfloating D/A converters or Digital outputs are hooked to the bus. Each D/A converter and each box with 16 digital outputs needs a proprietary adress. The oldest generation of boxes decodes only 5 adress bits. This has to be changed as soon as new boxes are hooked up, since the adressspace is used up.

```
unit NDYAG;
```

```
uses Serial;
```

Basic input/output for Nd:YAG laser using serial port.

```
unit Params;
```

```
uses Register,UniMess,RFSweep,crt,strings,tingu,exp,ShutDown;
```

This unit contains all Parameters used by the experiment and registers them, so that they can easily be stored, loaded, associated with the automatically created parameter menus. For each parameter an unique identifier for reference in \*.dat parameter files, a parameter range, a menu question and the color of the menu line have to be declared. New menu pages can be declared, they need a title and a link to a procedure that is

called after the menu is quit to eventually update settings of the experiment (used eg. for MOT frequencies). These procedures are declared in EXP.PAS.

Some nonregistered constants are also declared, like the time to wait after an IGBT command, maximum currents for the different power supplies. Space is reserved for storage of momentaneous output currents and tensions.

If a boolean type variable is registerd with the color "MarkColor", In the mainmenu a green spot appears at the side of the related parameters menu, if the option is switched on. This gives a overview at a glance of the momentary used experimental sequence.

```
unit PortList;  
uses crt,timing,io,MultiOut;
```

Sometimes during the experimantal sequence the ~30s calculation effort after a IOUpdate command to figure out which ports have to be updated is too high. Then a Portlist can be used. A Portlist is constructed by executing the sequence without actually writing to the ports. All I/O operations and WAIT commands are stored in the list. Then the list is written out "without thinking", which is faster.

A Portlist is usually started with the

```
SequenceStartPortlistDuringWait(10);
```

command, where 10 stands for 10ms preparation time and can eventually be longer. The operations are finally written when the

```
SequenceStopPortlist;
```

command is met. A Portlist should not be used for heavy traffic on the bus, like sweeps of magnetic fields, but for small short sequences like optical pumping, hyperfine transfer or image taking which are timecritical.

```
unit PulsGen;  
uses IEEE,crt;
```

Basic commands for the quantum composers 9300 series puls generator using IEEE.



```

program readHP8563A;
uses HP8563A;
    Main program for readout of HP8563A. Not actual part of LITHIUM.

unit Register;
uses Menu,serial,strings,crt;
    This unit contains routines to manage (store, load, send by serial
    port) the parameters of the experiment which are used by Param.pas
    which contains the actual declarations of the variables used in
    LITHIUM.

unit RFSweep;
uses crt,graph,Sequence,ExpIO,SigGen,params,strings,Register;
    class to register in parameterlist, modify graphically or with a
    numeric menu and execute an RFSweep.

PROGRAM SaveLeCroy;
USES CRT,LeCroy;
    Program to save trace from LeCroy, has nothing to do with LITHIUM

unit sequence;
uses crt,timing,portlist;
    This unit is used to do the reordering of the commands which are
    transmitted to the output ports. This is useful for example if a
    shutter has a delay time. Actually this reordering is never used
    in the momentary experimental sequence since it can not be used
    with the SequenceBlock or Portlists and thus is difficult to use.
    What is really used and useful to day is the following: The
    experimental sequence is executed twice. The first time nothing is
    written on the output ports, but general informations about the
    experimental sequence are gathered, like the total time for the
    sequence and the path followed by the atoms. Also the
    pulsgenerator is programmed and all its configurations used during
    the sequence are stored. (This is useful since it is slow in
    communication, but it can recall its parameters fast during the

```

actual execution.) The total time is needed in Vision to put the Andor camera in the external trigger mode just only about 100ms before the trigger arrives. (This is necessary because in this mode the CCD chip isn't well cleaned anymore. And during the sequence communication with Vision would be too slow.) The first execution takes about half a second. The second execution really writes on the ports.

During the first execution the experimental sequence is stored in memory, so reordering can be achieved. Since memory isn't sufficient to store every single command, in MagTrap.pas nearly every experimental step is stored as a block and only the block is stored. In fact, one could also completely switch off this memory storage option.

```
unit Serial;  
uses debug;
```

The serial communications unit in use in LITHIUM. Two units are at disposition for this. An old 9600 Baud version (serial1.pas) which is used right now and a modern 115200 baud version (serial2.pas) which I tried once but which produces too many missing characters during the transmission. An error correction protocol between LITHIUM and VISION would have to be implemented to use the latter. But since the communication load is not too high, the 9600 baud don't slow down the experiment but once when initially all parameters are transferred. After the initial synchronisation transfer only the few modified parameters are transferred from run to run.

```
unit Serial1;  
  see Serial
```

```
unit Serial2;  
  see Serial
```

```
unit Setups;
```

uses crt,menu,params;

in the SETUPS.DAT file different setups can be manually configured. Each setup consists of a subset of parameters and their new values when this setup is loaded. An usefull example is the calibration for the different lens systems. After change of the lenssystem the corresponding setup is loaded. To change a setup SETUP.DAT has to be modified. In contrast to Load and Save configuration the advantage is that only a subset of parameters is modified.

unit ShutDown;

uses Params;

After successull storage of the Parameters on harddisk a flag file is created. After loading the Parameters this file is destroyed. When the program terminates whith a crash the nonexistence of this file signals the problem during the next start of the program. At that point VISION and LITHIUM are probably out of sync. In this case a complete parameter retransmission to Vision is forced to resynchronise.

unit SigGen;

uses params,crt,IEEE,sequence,ExpIO;

This unit is used to program all synthesizers of whatever type in the same way. Like this it is easy to exchange one synthesizer by another, just change the type and GPIB number in the declaration of this file.

unit Strings;

To save memory space strings for the Parameter menus are stored on the heap in spaces exactly the size of the constant string. All the same each string is unfortunately stored twice since in pascal the adress of a constant string declared in the code can't be accessed.

unit Timing;

uses crt,dos,debug;

This Unit implements a good delay routine using the DAS 1601 Counter 0. If the DAS1601 card isn't detected the PC internal timer is used. The computer is synchronized to the counter to increase the reliability. startCounter has to be called shortly (less than the first wait-time and less than half a second) before the first Wait command. To identify Wait commands in the sequence to each delay time a femtosecond delay is added which of course doesn't play a role. But the amount of femtoseconds identifies the wait.

```
unit unimes;
```

```
uses Register,menu,serial,crt,strings,Exp,Vision;
```

This unit contains all the routines to manage universal measurements. To display the corresponding menus, store and load the data to and from UNIMESS.DAT. The routines necessary to actually execute the measurement are in Exp.Pas.

```
unit Vision;
```

```
uses Unimes,params,crt,serial,register,debug,ShutDown;
```

This unit contains the routines used to communicate with VISION. After the first experimental sequence is executed the first time with blocked outputs, some initial parameters like, camera, imagemodus, traptyp, the way of the atoms and the total time of the sequence are transmitted. Vision will wait till about 100ms before the end of the sequence and then put the camera into the external trigger modus. After the image has been taken, the parameters which changed since the last run are transmitted. If the COM port is changed it has to be modified here.